

Order form for Streptamer Custom Service

Last date of revision
June 2008



Version DC07-0001

Please fill out the following and fax to **(Europe) +49 (0) 551-50672-181**
(US) 1-888-531-6813

Streptamer Custom Service order	Please indicate
Production of a MHC I-Strep of the following species	<input type="checkbox"/> human <input type="checkbox"/> mouse (Cat. No. 6-9010-000)
MHC I heavy chain allele? (e.g. HLA-A*0201 or H-2 K ^b)	_____
Designation of antigen? (e.g. CMV PP65)	_____
Antigen peptide sequence? (8-11 amino acids in one-letter-code)	_____
Who will supply of peptide?	<input type="checkbox"/> I will supply the peptide (50 mg, HPLC purified, >80% purity, lyophilized, certificate of analysis will be provided, please refer to shipment instructions below) <input type="checkbox"/> IBA shall synthesize the peptide (Cat. No. 6-9020-000)
Has this combination of MHC I and peptide been successfully produced already?	<input type="checkbox"/> no, not to my knowledge <input type="checkbox"/> yes, in our lab <input type="checkbox"/> yes, it is offered elsewhere <input type="checkbox"/> yes, please refer to this publication _____
Initial order of custom made MHC I-Strep: (Note that for initial orders a discount of up to 25% is allowed)	<input type="checkbox"/> 50 tests; Cat. No. 6-9030-005 <input type="checkbox"/> 150 tests; Cat. No. 6-9030-015 <input type="checkbox"/> 500 tests; Cat. No. 6-9030-050 <input type="checkbox"/> 1000 tests; Cat. No. 6-9030-100
Order of the following Streptamer staining products: (Note that for initial orders a discount of up to 40% is allowed)	<input type="checkbox"/> Strept-Tactin PE 50 tests; Cat. No. 6-9040-005 <input type="checkbox"/> Strept-Tactin PE 150 tests; Cat. No. 6-9040-015 <input type="checkbox"/> Strept-Tactin PE 500 tests; Cat. No. 6-9040-050 <input type="checkbox"/> Strept-Tactin PE 1000 tests; Cat. No. 6-9040-100 <input type="checkbox"/> Strept-Tactin APC 50 tests; Cat. No. 6-9045-005 <input type="checkbox"/> Strept-Tactin APC 150 tests; Cat. No. 6-9045-015 <input type="checkbox"/> Strept-Tactin APC 500 tests; Cat. No. 6-9045-050 <input type="checkbox"/> Strept-Tactin APC 1000 tests; Cat. No. 6-9045-100
Order of the following Streptamer magnet assisted cell sorting products: (Note that for initial orders a discount of up to 40% is allowed)	<input type="checkbox"/> Streptamer Magnetic beads 5 preps; Cat. No. 6-9050-005 <input type="checkbox"/> Streptamer Magnetic beads 25 preps; Cat. No. 6-9050-025 <input type="checkbox"/> Streptamer Magnetic beads 75 preps; Cat. No. 6-9050-075 <input type="checkbox"/> Streptamer Magnetic beads 250 preps; Cat. No. 6-9050-250

Customer contact information		
	Billing Address	Shipping Address (if different from billing address)
Institution		
Department		
Address		
City		
Zip code		
State		
Country		
Contact First Name		
Last Name		
Email		
Phone		
Fax		
P.O.#		
Institution Type	<input type="checkbox"/> Academic <input type="checkbox"/> Industrial <input type="checkbox"/> Government	
How do you like to receive information, e.g. progress reports, results and documentation?	<input type="checkbox"/> by email <input type="checkbox"/> by fax <input type="checkbox"/> by normal post mailing	

Shipment of your peptide:
Please announce shipment of your peptide by email to service@iba-go.com prior to shipping. Please indicate the P.O.# on your sample shipment.

<ul style="list-style-type: none"> • Please inquire for exact delivery time. • "Tests" refers to staining with fluorescently labeled Strep-Tactin (1 test is sufficient to stain 5×10^6 T-cells) while "preps" refers to T-cell purification with Streptamer • Magnetic Beads (1 prep is sufficient to purify T-cells out of a total of 2×10^7 cells). • Please note that MHC I-Strep products are not available in North America. • IBA shall be entitled to take the produced MHC I-Strep into stock 6 months after delivery. • In case the ordered MHC can not be produced with the requested peptide the service fee plus eventually including peptide synthesis fee will be charged though. 	
Date: _____	Signature: _____